Chirurgische Standard- Instrumente / Surgical Standard Instruments

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Operating Manual
Page 18

DEUTSCH

ENGLISH

Inhalt

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Contents

1	Symbols	18
2	Introduction	18
3	Scope	18
3.1	Inteded Use	18
3.2	Contraindications.....	23
4	Warning.....	23
5	Handling	23
6	Preparation.....	24
6.1	Reutilization restrictions	24
6.2	Information on instrument preparation	24
6.3	Preparation at the place of use	24
6.4	Ultrasound bath (optional)	24
6.5	Manual cleaning	24
6.6	Mechanical cleaning.....	25
6.7	Drying	27
7	Sterilization.....	27
8	Maintenance, inspection	27
9	Packaging	27
10	Storage.....	27
11	Warranty / Repair	28
12	Manufacturer and servicing address.....	28

1 Symbols

Symbol	Definition
	CE marking
	Attention
	Manufacturer

2 Introduction

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following. In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

3 Scope

The use of Standard Surgical Instruments is mainly in the following areas:

- Dentistry
- Otolaryngology
- Gynecology
- Surgery

3.1 Inteded Use

Curettes:

curette; Curettes, titanium; Lupus curettes, -spoons, -scrapers: A spoon-shaped surgical instrument for use in gastroenterological and urological surgical procedures. This instrument group includes e.g. Biopsy spoon, ovule spoon, placenta spoon, uterine spoon, secretion spoon. It is a reusable instrument.

Biopsy curettes, Uterine curettes: A heavy, hand-held, surgical instrument with an elliptically windowed end on a long shaft for scraping out the uterine lining, scraping tissue during a gynecological procedure, and removing waste from the uterus. The instrument is usually made of stainless steel and can be either rigid or flexible. The working end is prepared sharp or blunt in different versions. The general product group includes uterine, placental and fallopian tube curettes. It is a reusable instrument.

Ear curettes; Paranasal Sinus Curettes: A surgical hand instrument having an annular, fenestrated, or spoon-shaped working portion which may be either sharp or blunt to provide tissue by scraping or scraping during an ENT procedure. It is a curette for general use in ENT surgery, it is made of stainless steel and is available in various designs and shapes. It is a reusable product.

Ring curettes: A surgical hand tool with an annular windowed working part for use in orthopedic surgery. It is used to remove bone tissue. It is a reusable product.

Exenteration scoops; Chalazion Spoon: An ophthalmic surgical instrument with a windowed, spoon-shaped or ring-shaped end that is either sharp or blunt. It is used to scrape out eye tissue. It is usually made of stainless steel. It is a reusable instrument.

Elevators:

elevators; Elevators: A surgical hand instrument for lifting, positioning or levering of bone structures, other anatomical structures or surgical material during an orthopedic procedure. It is made of stainless steel and is available in a variety of sizes, styles and working ends. It is usually a sturdy instrument with a blunt working end that can be pointed, rounded, flat, curved or hooked. The handle at the proximal end is also available in different versions. Some models have special functions, e.g. for ankle surgery, or to treat a dislocation. It is a reusable instrument.

Septum and Periostic Elevators: A surgical hand instrument for lifting, positioning or rupturing anatomical structures or surgical material in the ENT area. It can be a single or double instrument with a blunt end (leaf or hook in different versions) and a handle in the middle. Size and design of intended use, e.g. a dull impeller elevator, a septum elevator, a stronger nose artery elevator, a special viable ear elevator. The instrument is usually made of stainless steel. It is a reusable instrument.

Elevatorium: A surgical hand instrument for levering out, lifting and positioning a bone plate during surgery on the skull. Also known as bone flaps Elevator there is typically made of one piece and is robust, long, rounded, with a blunt blade at the distal working end and a handle at the proximal working end. It is typically made of high quality stainless steel. The instrument is reusable.

Tissue dissectors: A surgical hand instrument, usually in the form of spoon or rounded at the working end. It is used to separate soft tissue or body structures in interventions of general or plastic surgery. It usually has a handle that continues into a shaft that has a tip at the distal end. The tip may be pointed or flat, sharp or dull, angled or straight on the shaft. The instrument is made of high quality stainless steel and available in various shapes and dimensions. The instrument is reusable.

Comedo Extractors: A surgical hand instrument for levering out, lifting and positioning anatomical structures or surgical materials during ear, nose and throat (ENT) procedures. It may be a single- or double-ended instrument with blunt working end or hooks in various shapes without a cutting edge at the distal working end (s), and a handle positioned centrally or at the proximal end. Its shape and size depends on the intended use, for example: a blunt tympanic elevator, a rhinological or septal elevator, a strong nasal fracture elevator, or a specially shaped ears elevator. It is usually made of stainless steel. It is a reusable instrument.

Eye Specula: A surgical ophthalmic instrument for holding the edges and associated tissue of a cut in the eye during a surgical procedure. It usually has a self-retaining arm with curved hooks or leaves at the distal end that grip and hold the tissue. The proximal end is usually connected and has a lockable mechanism to hold the working end in an expanded position. It is made of stainless steel and is available in different sizes and designs. It is a reusable instrument.

Uterine depressors: A hand-held surgical instrument for holding off the uterus in order to facilitate an examination of the surrounding organs / tissues. It is typically a double-ended product in various designs (e.g., curved, spoon-shaped, or fenestrated) with a slightly angled working end to better manipulate the organs. It is made of stainless steel and is a reusable instrument.

Hooks: (□hauthäkchen)

Tissue Hooks; Micro Hooks: A shaft-like surgical instrument, having a proximal end forming a handle which can be formed in various embodiments. The shaft at the distal end merges into one or more parallel hooks that apply tension to the skin. It is a reusable instrument.

Hooks; Exploration-& Coagulation Hooks; Root retractors; excavators; excavator- curette: A shaft-like surgical instrument, having a proximal end forming a handle which can be formed in various embodiments. The shaft at the distal end merges into one or more parallel hooks to hold tissue. It is a reusable instrument.

Trachea retractors; Thyroid Gland retractors; Hooks: A surgical instrument for stably holding the trachea during a tracheostomy and / or separating / pushing away the margins of a tracheostomy to set a tracheostomy tube. The instrument may be of various types, e.g. Single or double ended, with a wide retractor on one side and opposite a pointed hook or several parallel hooks, which can be bent at different angles. It is made of stainless steel. It is a reusable instrument.

Hooks for posterior nares: A shaft-like surgical instrument whose proximal end forms a handle that can vary in shape. The shaft tapers toward the distal end toward a hook or hooks used in plastic nose surgery. In the multi-hook variant, the hooks are arranged in parallel. The instrument is reusable.

Ear Hooks; paracentesis Needles; Foot plate instruments: A fine surgical instrument that consists of a handle that tapers into a shaft and ends in a fine, 90 ° bent tip. It is used during ENT surgical procedures to repair structures in the middle ear, e.g. to edit the stapes. It is a reusable instrument.

Eye Hooks; Strabismus Hooks: An ophthalmic surgical instrument with a shaft-like handle in various versions that tapers towards the working end. The distal part of the instrument is blunt or pointed bent or angled to better manipulate the ocular tissue or to remove debris. It is a reusable instrument.

Laser Instruments; Hooks: An instrument used in conjunction with a suitable endoscope during endotherapy. It is used for mechanical work, e.g. grasp tissue or foreign objects, to bring medical fluid or width of an area / lumen. It works without electricity, including e.g. High-frequency, electromagnetic, ultrasonic or laser energy. It is a non-active product that is reusable.

Sharp Hooks: A shaft-like surgical instrument whose distal end is in the form of a hook or more parallel hooks for fixing or pulling uterine fibroids during a gynecological procedure. It is typically a one-piece stainless steel instrument with a handle in various designs at the proximal end. It is a reusable instrument.

Dilators

Uterine-Dilators: A round rod-like massive surgical instrument that is inserted through the cervix and dilates the cervical canal. Commonly known as a uterine dilator, it is typically made of metal or plastic and has a tapered rounded distal tip. It can have one or two different diameters, one at each end of the instrument, and can be part of a set of size-increasing diameters which are inserted successively. The instrument allows insertion of a curette or other instrument into the uterus and / or expands the uterus prior to surgery by pushing against the walls of the uterus from the inside. The instrument is reusable.

Trachea Dilators: A hand surgical instrument used during tracheal surgical procedures to dilate tracheal structures e.g. during a tracheostomy and / or to assist in the insertion of a tracheal cannula. It can have different designs and sizes, e.g. a slender hollow or solid body made of metal, plastic or other suitable materials in cylindrical form, or a scissors-like instrument that activates by closing the two ring handles two or three blunt branches in opposite directions. The instrument is reusable.

Bile Duct Dilators: A surgical instrument consisting of a slender hollow or solid body of metal or plastic in a cylindrical shape in a variety of sizes for dilating the bile duct, e.g. at the union of the ductus sympathicus with the ductus cysticus. The instrument is reusable.

Prostatic Dilators: A round, rod-like, surgical instrument for examining and dilating the urethra. It has a semi-rigid, conical shaft that merges into a narrow, bulbous finish. The instrument can be made hollow to be able to pass other, narrow instruments. It is a reusable instrument.

Vascular Dilators: A surgical instrument that is inserted into a vessel and used to dilate or adjust a vessel during a surgical procedure. The instrument is not intended for use on the blood vessels arteriae pulmonales, aorta ascendens, arcus aortae, descending aorta to bifurcation aortae, coronary arteries, common carotid artery, external carotid artery, internal carotid artery, cerebral arteries, brachiocephalic trunk, venae cordis, vena pulmonary, superior vena cava, and inferior vena cava. It is a rod-like instrument with a long, slender solid body and a distal tip in various shapes and sizes. It is available in a variety of sizes and degrees of softness and is typically made of metal. The instrument is reusable.

Lachrymal Dilators: An ophthalmological instrument for use during ophthalmological procedures to dilate a cavity, canal or opening. Eye dilators also include instruments used for lacrimal or lacrimal canal surgery. It is a reusable instrument.

Needles:

Intestinal Needles; Muscular- & Skin suture; Reverdin Needles; Ligature Needles; post mortem needle: A cylindrically shaped, solid instrument of various diameters, whose distal tip can be made in various geometric shapes. It is used for introducing and pulling suture through superficial tissue, e.g. Skin or muscle, used. The suture serves to piece two wound edges together. It is a reusable instrument.

Eye Needles: A cylindrical, solid instrument of various diameters, whose distal tip can be made in various geometric shapes. It is used to insert and pull suture through ocular tissue. The seam connects two tissue surfaces. It is a reusable instrument.

Guide needles: A product for connecting a catheter to another product (e.g., a plug) or for introducing a flushing or drainage catheter into the body. It is a reusable product.

Spatula:

Tissue spatula: A surgical instrument, usually made of stainless steel, that is used to manipulate tissue or to remove material from a surface. The instrument consists of a proximal handle and a flat blade with no sharp edges at the distal end. The shaft or the blade of the instrument can be straight, curved or variously angled. It is a reusable instrument.

Plaster Spatula; Abdominal-& Intestinal Spatulas; Organ spatulas: A surgical hand instrument for manipulating the intestine and abdominal organs during an abdominal surgical procedure. It typically has a large, flat, blunt blade (without sharp edges) at the distal end and a strong handle at the proximal end. It is usually made of stainless steel and is available in various designs, e.g. curved and flat, or straight and flat with a rounded end, or right angled. Some models have a ridged leaf surface to better hold the slippery abdominal tissue. The handle is typically designed so that it can be kept non-slip. It is a reusable product.

Lung spatulas: A surgical instrument, usually made of stainless steel, which is used to manipulate lung tissue surfaces. The instrument consists of a proximal handle and a flat blade with no sharp edges at the distal end. The shaft or the blade of the instrument can be straight, curved or variously angled. It is a reusable instrument.

Blades: A surgical hand instrument, usually in spoon shape or rounded at the end of work. It is used to separate soft tissue or body structures in interventions of general or plastic surgery. It usually has a handle that continues into a shaft that has a tip at the distal end. The tip may be pointed or flat, sharp or dull, angled or straight on the shaft. The instrument is made of high quality stainless steel and available in various shapes and dimensions. The instrument is reusable.

Bladder Retractors; Gall Stone Scoops; Uterine Curettes: A spoon-shaped surgical instrument for use in gastroenterological and urological surgical procedures. This instrument group includes e.g. Biopsy spoon, ovule spoon, placenta spoon, uterine spoon. It is a reusable instrument.

Tongue Depressors: A surgical hand-held tongue removal instrument to facilitate examination of the surrounding organs / tissues. It is typically a double-ended, flat product of various sizes with a slightly angled working end to better manipulate the tongue. It is made of stainless steel and is a reusable instrument.

(□Zementspatel)

A hand surgical instrument for applying orthopedic bone cement to a surface. It is usually made of stainless steel and consists of a proximal handle and a flat blunt blade (no sharp edges) at the distal end. The shaft or the blade of the instrument can be straight, curved or variously angled. The instrument is reusable.

Specula:

Ear Specula: A surgical instrument, usually composed of two articulated branches, available in a variety of sizes, shapes and contours. It has a blunt distal end and is used in insertion and opening to stretch or stretch tissue in the opening of the ear. It is used for examination or during procedures. It is usually made of stainless steel with a self-opening mechanism or an adjustable mechanism. It is a reusable instrument.

Rectal Specula: A surgical instrument, usually composed of two articulated branches, available in a variety of sizes, shapes and contours. It has a blunt distal end and is used for insertion and opening to stretch or stretch the rectal opening (canal). It is used to study or introduce other instruments. It is usually made of stainless steel. The portion of the branches in the proximal area forms the handle while the distal area opens when you squeeze the handle. It is a reusable instrument.

Mirrors:

Laryngeal Mirrors; Laryngeal mirror handles; Middle Ear Mirror: A surgical hand instrument with a polished surface, for reflecting sufficient undiffused light to obtain a visible image of an object during an ENT procedure. This mirror is mounted on a long, narrow handle, which is held by the doctor in the vicinity of the object to be observed. It is a reusable instrument.

Others:

Mallets: A manual surgical instrument for applying force (e.g., by beating) to another product (e.g., a chisel, pry bar, driver) during a surgical procedure. It is available in different designs and materials. The shaft is usually made of metal (chromed iron or stainless steel), but may also be made of synthetic materials (e.g., tufnol). The hammer head is usually solid, annular, mostly double-ended,

and is made of solid (iron / steel) or a shock-absorbing material (e.g., rubber, plastic, or tufnol). It is a reusable instrument.

Myoma Screw: A manual surgical hand instrument for detecting and holding a fibroid (i.e., a tumor of muscle tissue, usually in the uterus) during its surgical removal. It is typically a one-piece instrument, similar to a corkscrew with a pointed, spindled working end (for screwing into the myoma) and a long, slender shaft and a T-shaped handle. It is usually made of stainless steel. It is a reusable product.

Cotton Applicators: A narrow, rod-like, flexible metal surgical instrument, usually with a blunt, beaded, or pointed tip, for examining fistulas, lumens, or wounds during an ENT procedure. The distal end, in front of the bead / tip, may be bent at right angles to the shaft. It is a reusable product.

Percussion Hammers: A surgical hand instrument for gently tapping the reflex points (e.g., near the knee, ankle). It is composed of a handle and a shaft, usually made of stainless steel. Its head may be annular, wedge-shaped or dome-shaped and made of a soft material (e.g., rubber or plastic). Some designs may be provided with a spike, brush or other interchangeable components for neurological examinations. It is a reusable instrument.

(□Dilatationsbougies)

A round, rod-like, surgical instrument for examining and dilating the urethra. It has a semi-rigid, conical shaft that merges into a narrow, bulbous finish. The instrument can be made hollow to pass other, narrow instrument can. It is a reusable instrument.

Table top cutters: A heavy duty surgical hand instrument with specially designed powerful handles and a cutting jaw for cutting orthopedic wire, small diameter pins or cerclages. It has a scissor-type design with straight or curved handles and is made of stainless steel. It is available in various sizes and the cutting edges are usually short, strong and dull. The final part consists either of a single or double ratio screw connection to apply a greater force. The instrument, commonly referred to as a wire cutter, may be designed as a front or side cutter, some models have carbide inserts to extend shelf life. It is a reusable product.

Tuning forks: A U-shaped instrument, usually made of stainless steel, with a handle at the base of the "U". The vertical forks of the "U" are cut to a specific length to produce a sound at a particular wavelength when the "U" is struck against a hard object, usually rubber. It is typically used during an ENT examination to test the auditory sharpness. It is a reusable instrument.

(□Areola-Marker)

A hand-operated surgical instrument used for removing blemishes or for sampling skin tissue. It usually consists of a handle with tips at the end of work in various configurations; some peaks include a closed cutting zone, edged by a cutting blade. This is a reusable instrument.

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A surgical instrument with a sharp pyramidal or conical tip used to pierce body cavity. It is usually assembled and used with a compatible sleeve, filling its lumen, thereby allowing insertion of the instrument. After piercing, the trocar is withdrawn and provides a working channel into the body cavity. It is a reusable instrument.

Screwdriver: A tool, which can also be called a surgical instrument, with a shaft that fits into a screw head to tighten or remove this screw. The working end may e.g. for a slotted screw, Phillips screw, Phillips, Hexagon. The shaft end may have a handle for manual operation or a connection for a chuck for electric drive. It can have an indication of the torque. It is a reusable instrument.

Lead Hands: A personal protective equipment used to shield the entire hand of an operator or other from unnecessary primary and scattered radiation in the context of diagnostic and therapeutic procedures. Each finger is shielded individually. When used for X-ray and nuclear medicine, they usually consist of a liquid resistant outer covering surrounding a thin layer of lead or lead equivalent material. The gloves provide your fingers with radiation protection from all directions, the top and bottom of the hand. Some offer an additional physical barrier to protect against liquid radioactive materials and body fluids. It is a reusable instrument.

Extension Bows: The product is part of a traction unit. It consists of metal bars with a curved frame, for attaching hooks or pins for exerting a pull on the patient. The traction unit exerts a tensile force to allow the body parts, e.g. the cervical or lumbar spine, by applying a strap around the head or pelvis. The adjustable spreader bar adjusts the angle at which the belt exerts the tension. It is a reusable instrument. **Endarterektomiestripper; Varicose vein cutlery; Venous Instruments, Venous Hooks, Venous Extractors:** A surgical instrument to cut a piece out of a vessel. A stripper can be used to excise part of a vein or artery except for the arterial arteries pulmonary, ascending aorta, aortic arc,

descending aorta, aortic bifurcation, coronary artery, common carotid artery, external carotid artery, internal carotid artery, arteriae cerebrales, brachiocephalic trunk, venae cordis, vena pulmonales, superior vena cava, and inferior vena cava. Vessel Stripper is available in two different versions: 1) a flexible stainless steel cable with a peeling dome or disc at one end and a guide tip at the other end; and 2) a rigid rod terminating in a closed loop or loop, e.g., an external stripper. The instrument is reusable after proper cleaning has been done.

Measuring probes: A product / instrument used for comparative measurement in clinical use, e.g. inner and outer diameters, lengths, depths or thicknesses to measure. It is a non-calibrated product that is reusable.

Ligature Guides; Wire strut guide: A surgical instrument for guiding wire or ligatures through tissue. Typical guides may have flexible or rigid working parts with hooks, clips or claws for holding the material while passing through the tissue. The working end can have a button or an eye.

Füllungsinstrument: A dental instrument for applying filling and restorative material to the patient's mouth. It is a reusable product.

Sharp spoon: A spoon-shaped surgical instrument for use in gastroenterological and urological surgical procedures. This instrument group includes e.g. Biopsy spoon, ovule spoon, placenta spoon, uterine spoon, secretion spoon.

3.2 Contraindications

There are no known contraindications.

4 Warning

!	Medical products are delivered in a non-sterile condition and must be cleaned, disinfected and sterilized prior to their initial use.
!	The use of faulty instruments is in principle forbidden and they have to go through the whole cleaning process before return.
!	Please take into consideration that through higher power a bigger damage of the tissue can result: f.e. on forceps: the power at the end of the jaw is higher than at the tip of the jaw
!	Please observe the additional information enclosed with the products.
!	Remove all protective sleeves and films prior to first using or preparation for use.
!	The safe combination of different products or of products with implants must be reviewed prior to clinical application by the user.
!	Avoid improper throwing or dropping of instruments
!	To prevent all contact corrosion, instruments with damaged surfaces must be separated immediately.
!	If used on patients with Creutzfeldt-Jakob disease or HIV infection, we are unable to accept any responsibility for reuse.

5 Handling

The type of treatment must be determined in each individual case by the surgeon in cooperation with the internist and the anaesthetist.

For operational use in various surgical disciplines must be done by appropriately trained and qualified personnel.

6 Preparation

The person in charge of preparatory treatment is responsible for ensuring that the treatment is duly carried out using the relevant equipment, materials and personnel in the treatment facility and so achieves the desired result. This necessitates validation and routine monitoring of the process used. We urge you to take note of the national regulations dealing with instrument preparation.

6.1 Reutilization restrictions

Frequent repeat preparatory treatment has minimal effects on the product. The end of the product life is normally determined by wear and damage due to use

6.2 Information on instrument preparation

- Use cleaning and/or disinfection agents with a pH-value within 9-10. Please observe manufacturer instructions regarding dosage, exposure time and renewal of solutions.
- Do not use hard brushes or coarse abrasive cleaners.
- Never leave instruments in cleaning or disinfection agents for longer than the specified time.
- Only used demineralized water for rinsing.
- Rinse and dry carefully through channels and pipes.
- Sensitive instruments must be cleaned in a storage or clamping fixture.
- Observe manufacturer instructions of cleaning – and sterilizing equipment.

6.3 Preparation at the place of use

Directly after using remove coarse dirt of the instruments and rinse out the working cannulas. Do not use fixing agents or hot water ($> 40^{\circ}\text{C}$), as this results in residues becoming fixed and can affect the success of the subsequent cleaning operation

Dismantle and/or open instruments as far as possible. Within short time after use the instruments clean the instruments for reducing a drying of the residues. **Use cleaning and/or disinfection agents with a neutral pH-value <10 .**

This enables an easier cleaning. If instruments come into contact with corroding medicines or cleaning agents, wash these up with water immediately after use.

6.4 Ultrasound bath (optional)

All instruments must be opened, dismantled and any cavities rinsed through.

Place instruments in the screen basket in such a way that overlaps and contact between instruments are avoided. Add cleaning

agent to the water and adjust the temperature of the solution in line with the cleaning agent manufacturer's instructions.

The cleaning in the ultrasound bath should be at 35-40 kHz, 5 minutes at least.

Subsequently rinse instruments including all cavities before cleaning and disinfection.

Medicine products which are possessing a bad ultrasound transmission, e.g. soft material are not usable for the ultrasound bath.

6.5 Manual cleaning

We recommend using lint-free soft cloths, plastic brushes and/or cleaning guns for manual cleaning. The cleaning agent should not have any fixing properties. The cleaning fluid must be changed daily.

Hints:

- If possible, apply the cleaning agent into cavities at pressure and carefully brush out.
- Be careful with the cleaning of barriers, slits, tight entries and joints etc.
- It has to be secured that instruments with joints will be cleaned in an opened and closed position.
- Use adequate brushes for the cleaning of cannulas, working channels and sack holes
- Flexible (complexe) instruments with unvisible surface and chuck have to be cleaned manually.
- We recommend for all push shafts, -and pipe shafts instruments and instruments whose surfaces are on top of each others during the cleaning (f.e.bone forceps and gouge forceps) a manually cleaning for an optimal cleaning result without residues

After that rinse the instruments with fully demineralized water extensively. Subsequently dip the instrument in disinfectant and use suitable cleaning instruments such as synthetic soft brushes. Observe the cleaning agent manufacturer's instructions. At the end of the cleaning process, repeat the rinsing action with fully demineralized water and dry the instruments completely.

6.6 Mechanical cleaning



On the basis of international standards (EN ISO 15883) and national directives, only validated machine cleaning and disinfection methods may be used. For the mechanical cleaning we recommend a standard programme for surgical instruments, f.e. instruments from Miele.

Observe the following by loading:



- Place the dismantled/opened instruments securely in the tray.
- Instruments with openings and gaps have to be faced down with the opened side so that they can be cleaned and no water of the cleaning process is collecting inside them. If available use balanced devices for rinsing

- Place the instruments with joints in an opened position into the cleaning,-and disinfection machine.
- If needed use an adapter for the cleaning

- Do not overload trays, avoid creating any overlaps.

Preliminary rinsing (cold, if applicable fully demineralized water without additives) is followed by chemical.

The chemical cleaning should take place at 40°C -60°C for at least 5 minutes.

We recommend products with a **pH-value within 9-10**, e.g. MediClean forte from Dr.Weigert. The cleaning agents used should be selected depending on the material and properties of the instruments and in accordance to national regulations: If there is a high chloride concentration in the water, pitting and tension crack corrosion can occur on the instruments. The occurrence of this type of corrosion is minimized by using alkaline cleaning agents and demineralized water. By adding an acid-based neutralization medium, the rinsing off of alkaline cleaning agent residues is facilitated during the first intermediate rinsing process (warm or cold water).

In order to prevent the formation of deposits, it is advisable to use neutral cleaners where the water quality is unfavourable. After the second intermediate rinsing process, thermal disinfection takes place.

The thermal disinfection should take place at temperatures of between 80 and 95°C, with an exposure time as outlined in EN ISO 15883.

After the finishing of the programme take the good out of the machine because corrosion can arise if the instrument remains in the machine.

6.7 Drying

Ensure adequate drying by the cleaning and disinfection device or using other suitable measures.

The drying temperature should not exceed 95°C to avoid material-related ageing processes.

7 Sterilization

Prior to sterilization, products must undergo cleaning and disinfection, be rinsed off without residue using demineralized water and subsequently dried. HEBUmedical recommends using a validated steam sterilization process (e.g. sterilizer in compliance with EN 285 and validated in accordance with DIN EN ISO 17665-1).

On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132° C)** with a **minimum dwell period of 4 minutes** and subsequent vacuum drying.

The vapour must be free of ingredients, recommended limiting values of feed water and vapour condensate are determined through EN 285.

Other sterilization processes are compatible but not validated from HEBUmedical.

When loading, observe the recommended total weight. After the sterilization, check the sterile product packaging for damage, and inspect the sterilization indicators

8 Maintenance, inspection

After cooling to room temperature, the instruments must be visually inspected for protein residues and other contamination, paying particular attention to cavities, blocks, inclusions, pipes, and other inaccessible areas. Instruments which are not free of residues must be returned for a complete retreatment process.

Carry out the functional check mentioned above. Instruments with stains, which are blunt, bent, no longer function or which are otherwise damaged must be segregated!

9 Packaging

Compliant packaging of products for sterilization in line with ISO 11607. Packaging used must be suitable for the instruments and protect them from microbiological contamination during storage. The seal must not be under tension. HEBUmedical recommends container for sterilized packaging.

10 Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from 5°C to 40°C. Protect from the effects of the sun's rays and artificial light. Warranty / Repair

11 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

12 Manufacturer and servicing address

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany
Tel. +49 7461 94 71 - 0
Fax +49 7461 94 71 - 22
eMail: service@HEBUmedical.de
Web: www.HEBUmedical.de

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**Our vast range of products covers over 10 000 various instruments.
They can be found in Internet or request our catalog.**

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