Klemmen und Zangen / Clamps and Forceps

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Seite 4

Operating Manual

Page 17

DEUTSCH

ENGLISH

Inhalt

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Contents

1	Symbols	17
2	Introduction	17
3	Scope.....	17
3.1	Inteded Use	17
3.2	Contraindications.....	21
4	Warning.....	22
5	Handling.....	22
6	Preparation	22
6.1	Reutilization restrictions.....	22
6.2	Information on instrument preparation	22
6.3	Preparation at the place of use	23
6.4	Ultrasound bath (optional)	23
6.5	Manual cleaning	23
6.6	Mechanical cleaning	23
6.7	Drying.....	25
7	Sterilization	25
8	Maintenance, inspection.....	25
9	Packaging	26
10	Storage	26
11	Warranty / Repair	26
12	Manufacturer and servicing address.....	26

1 Symbols

Symbol	Definition
	CE marking
	Attention
	Manufacturer

2 Introduction

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following. In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

3 Scope

The use of forceps and clamps is mainly in the following areas:

- Thoracic Surgery
- Vascular Surgery
- tracheotomy
- Visceral Surgery
- Gynecology

3.1 Inteded Use

Bulldog Clamps; Bulldog Clamps, atrauma: A coarse toothed surgical instrument for grasping, squeezing, connecting or holding an organ, vessel or tissue. The instrument is usually made of stainless steel and is available in two versions: 1) a flat, tweezer-like design with legs that intersect between the jaw part and the gripping surfaces. When the user presses the gripping surface, the jaw part opens or closes. 2) A locking, scissor-like design with ring handles. Both versions are available with different mouth shapes: straight, bent or angled. It is a reusable instrument.

Bulldog Clamps, titanium: A coarse toothed surgical instrument for grasping, squeezing, connecting or holding an organ, vessel or tissue. The instrument is made of titanium and is available in two versions: 1) a flat, tweezer-like design with legs that intersect between the jaw part and the gripping surfaces. When the user presses the gripping surface, the jaw part opens or closes. 2) A locking, scissor-like design with ring handles. Both versions are available with different mouth shapes: straight, bent or angled. It is a reusable instrument.

Gall Duct Clamps: A surgical instrument for holding and manipulating the gallbladder during a surgical operation. It is typically a scissors-like, self-retaining instrument with ring handles. It is made of stainless steel in different sizes. The end of work is usually grooved. The branches are interlocked to find additional adhesion to the gallbladder. It is a reusable instrument.

Hysterectomy Forceps: A surgical instrument for grasping, pulling or pushing the uterus during a hysterectomy (removal of the uterus). It typically has a self-holding, scissor-type design with ring handles made of stainless steel. It is available in different sizes. The working end may be made in various designs, e.g. straight or bent. Some models may have an additional gripping tooth in their mouth. The mouth surfaces are roughly grooved, they can also be gerieft in the longitudinal direction in order to hold the organs better. The instrument is also known as a parametrium clamp or vaginal clamp. It is a reusable instrument.

Kidney Pedicle Clamps; Kidney and Pankreas Clamps: A surgical instrument for grasping and lifting the kidneys during a surgical procedure. It is usually designed scissor-like and self-locking, provided with ring handles. It is made of stainless steel and is available in different sizes. The end of the work has gripping industries shaped as very large, grooved, oval rings or semi-open circles. It is a reusable instrument.

Dissecting and Ligature Forceps: A surgical instrument for grasping, manipulating, squeezing or ligating tissue and ligatures during a dissection and / or autopsy. It is usually available in two versions: 1) a self-retaining, scissor-type design with ring handles and a working end with various designs, e.g. straight, curved or angled, with parallel grooves or teeth to improve grip; and 2) tweezer-like design (also micro / fine instrument) with two intermeshing branches with serrated tips at the working end. The instrument is made of stainless steel and available in different sizes. The working end can be reinforced with carbide inserts. The instrument is reusable.

Tubing Clamps: A surgical instrument for squeezing a tube during surgery or surgery to stop flow. It is typically a locking scissor-type instrument with stainless steel ring handles. It is available in various sizes with strong, thick, wide mouth parts, so as not to damage the hoses. Some models have scored mouth surfaces. It is a reusable instrument.

Sponge Holding Forceps; Cotton Swab Forceps: A scissors-shaped, surgical instrument with ring handles, whose working end is shaped as a ring, loop or ellipse to hold swabs or bandages. It is a reusable instrument.

Sinus Forceps; Clamp Forceps; Clamp Closing Forceps; Clamp holders: A scissors-type surgical instrument with ring handles. The leaves are serrated and either straight or curved. The serrated area has a central, non-toothed area in its length to prevent injury to the surrounding vessels. It is a reusable instrument.

Towel Forceps: A surgical instrument for holding wipes and other products, e.g. Cables, guides that must be securely attached near the surgical field, e.g. Towels covering the operating table. The instrument is typically made of stainless steel and may have different designs. It may, for example, have two branches crossed at the end with ring handles, or it may be a one-piece, A-shaped, end-connected instrument, the working end may be pointed, truncated, fenestrated or differently ridged. The instrument is commonly known as a cloth clamp or cloth clip. It is a reusable instrument.

Bronchus Clamps: A surgical instrument for atraumatic compression of the bronchi. It is a reusable instrument.

Blood Vessel Clamps; Sigmoid Hemostatic Forceps: A surgical device for achieving temporary hemostasis during surgery on vessels other than the blood vessels, arteriae pulmonales, ascending aorta, aortic arc, descending aorta, aortic bifurcation, coronary artery, common carotid artery, external carotid artery, internal carotid artery, arteriae cerebrales, truncus brachiocephalicus, venae cordis, vena pulmonales, superior vena cava, and inferior vena cava. It is a reusable instrument.

Approximators: A surgical device for the temporary retention of vessels other than the arterial vessels pulmonary, ascending aorta, aortic arcus, descending aorta, aortic bifurcation, coronary arteries, common carotid artery, external carotid artery, internal carotid artery, cerebral arteries, brachiocephalic trunk, venae cordis, vena pulmonales, superior vena cava and inferior vena cava. It is a reusable instrument.

Peritoneal Clamp Forceps: A surgical instrument for holding the peritoneum during a surgical procedure. It is a reusable instrument.

Multi-Purpose-Clamps; Polypus and Dressing Forceps; Nasal Polypus Forceps; Tissue Seizing Forceps with T/C jaws; Intestinal and Tissue Grasping Forceps; Capsule Forceps; Prostatic lobe forceps: A surgical instrument for atraumatically holding / gripping and / or squeezing intestinal tissue, tissue and some organs during a surgical procedure. It typically has two designs: 1) a self-holding, scissor-like design with ring handles; or 2) a large tweezer-like design with straight branches interlocking at the proximal end. It is available in different sizes. The working end may be made in various designs, e.g. ring-shaped, triangular, with parallel grooving or tothing to give grip, or inwardly bent profiles for grasping. The industries are typically wide and thin. It is made of stainless steel. It is a reusable instrument.

Intestinal and Stomach Clamp; Intestinal Clamp Forceps: A surgical instrument for atraumatically gripping, compressing, connecting or holding the intestine during gastrointestinal procedures. It is a reusable instrument.

Penis Seizing Forceps: A surgical instrument for grasping and holding the penis during a procedure. It is a reusable instrument.

Sterilizing forceps: A scissors-type surgical instrument with ring handles whose leaves at the end of the joint form open circles, loops or ellipses specially designed to grip and apply sterile instruments or implants directly from a sterilizer. This is a reusable instrument.

Haemorrhoidal grasping Forceps: A scissors-type, surgical instrument with ring handles, with triangular, toothed jaw part at the working end. For use in rectal / abdominal surgery. The instrument is reusable after sterilization.

Clip Applying Forceps: A surgical instrument adapted for the use of hemostatic clips for connecting blood vessels except than arteriae pulmonales, ascending aorta, aortic arc, descending aorta to aortic bifurcation, coronary artery, common carotid artery, external carotid artery, internal carotid artery, arteriae cerebrales, truncus brachiocephalicus, venae cordis, vena pulmonales, superior vena cava, and inferior vena cava. It is a reusable instrument. Forceps: A surgical hand instrument, consisting of two branches, which go over a conclusion in two jaws. The jaw parts are the actual working part of the forceps, which are used in various surgical, surgical treatments and are therefore available in different shapes. It is a reusable instrument.

Plaster Cast Instruments: A large, scissor-type, hand-held instrument whose leaves, distal to the end, merge into different mouth shapes suitable for cutting thick layers of plaster or synthetic material used for a plaster cast for various body parts. The tip of one of the distal blades of the scissor-type cutting mechanism is typically rounded at its lower edge to prevent injury to the patient when cutting the cast. The product is usually made of stainless steel and usually has a spring mechanism that helps to open the mouth. It is a reusable product.

Wire tightening-and twisting forceps; Wire tightener; wire and pin cutters; wire seizing forceps with tungsten carbide inserts; Wire cutting pliers: A surgical instrument with specially designed powerful jaws for holding, cocking and / or twisting wire introduced during a surgical procedure. It usually has a scissor-like design with ring handles, possibly with a lock and is made of stainless steel. It is available in different sizes and the working end can be made in different mouth shapes, e.g. typically short and wide with carbide inserts. Some may also be long, with side pins for attaching the wire. The jaw and trailer are usually heavily serrated. It is a reusable instrument.

Catheter Introducing Forceps: A hand-held instrument for gripping a tube (a catheter or an endotracheal tube) that is inserted or removed into the trachea. It also serves to grasp a foreign body from the trachea. It is commonly known as "Magill catheter insertion forceps" and has a scissor-like appearance with handholds. It is made of stainless steel. It is available in different finishes and the working end usually has gripping branches with narrow, round or S-shaped open mouth. The mouth surfaces are grooved to provide better grip. It is a reusable instrument.

Drill guide and protecting forceps: A surgical hand instrument in the form of a cannula for simultaneously placing and anchoring a rotating drill (i.e., a surgical twist drill) through hard tissue. Surrounding soft tissue is protected while the drill is in operation. It usually has ridges at the distal end to give a secure placement on the tissue surface. At the proximal end is a handle for holding and moving. The product is usually adjusted to the drill size and is made of stainless steel. It is a reusable product.

Flat nose pliers: A surgical instrument with specially designed sturdy handles and grasping mouth (usually parallel) for grasping and holding an object during a surgical procedure. It has a scissor-type design with curved handles and is made of stainless steel. It is available in different sizes and the jaw part is usually wide with a ridge. The branches are connected by a swivel joint or a double-pitch screw connection to allow greater force to be exerted. It is a reusable instrument.

Bending forceps: A surgical tool for bending wire that serves as an ossicular replacement in the ENT area. It is a reusable product.

Occlusive Clamps: A surgical instrument for atraumatically gripping, squeezing, connecting or holding the rectum or rectal canal. It is a reusable instrument.

Foreign Body Forceps: An instrument used in conjunction with a suitable endoscope during endotherapy. It is used for mechanical work, e.g. for grasping tissue or foreign, to introduce medical fluid or areas of a lumen. It works without electricity, including e.g. High-frequency, electromagnetic, ultrasonic or laser energy. It is reusable.

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Ethmoid Forceps: A surgical instrument for the removal of tissue or bone specimen during ear, nose and throat (ENT) procedures. The instrument can be equipped with a ring handle or pistol grip. The ring handles or pistol grips have extended shafts provided with punch-like or mechanical locking mechanisms at the distal end. The instrument is a reusable instrument.

Ovarian Elevating Forceps; Uterine Hemostasis Forceps: A surgical instrument used for general atraumatic grasping, pulling or compression of internal structures during a gynecological surgical procedure. It is usually designed like scissors with ring handles and is made of high quality stainless steel. It is available in various sizes, the working end may have a variety of leaf shapes, e.g. straight or bent. The leaves are usually flat in profile and thin, and rounded, so they do not hurt the internal organs. This is a reusable instrument.

Obstetrical Forceps: An obstetric surgical instrument to help with difficult births. It usually has two curved blades which are individually inserted and then assembled together like a handle which is placed around the child's head to pull or rotate to facilitate passage through the birth canal. It is made of stainless steel. It is a reusable instrument.

Septum Forceps; Nasal Septum Forceps; Goitre grasping forceps: A general surgical instrument for gripping, holding or manipulating anatomical structures during an ENT procedure at e.g. the bronchi, esophagus, trachea, larynx, pharynx, nose, or ear. It usually has two main designs: 1) a locking, scissor-type design with ring handles and a working end with different mouth designs, e.g. straight, angled, or curved with teeth or ridges to improve retention, and 2) a tweezer-like design (may also be a micro-fine instrument) with two connected legs with toothed jaw. It is made of stainless steel. It is a reusable product.

Vitrectomy grasping forceps: A surgical eye instrument for removing foreign bodies or the vitreous body from the eye. It is a reusable instrument.

Eye Forceps: A surgical ophthalmic instrument for holding, manipulating, compressing, pulling or connecting ocular tissue or surrounding tissue during a surgical procedure. It is a reusable instrument.

Tonsil Seizing Forceps: A surgical instrument for grasping, securing and manipulating the tonsils during an ENT procedure, usually during an almond removal. It typically has a scissor-type, self-retaining design with ring handles and is made of stainless steel. It is available in different sizes. The working end is typically gripping with curved blades of various designs and / or ending in oval rings or narrow serrated triangles. The leaves are usually coarsely toothed to better grip the tonsillar tissue. It is a reusable instrument.

Tongue Holding Forceps: A surgical instrument for grasping, holding or manipulating the tongue during a surgical procedure. It usually has a self-retaining, scissor-like design with ring handles. The working end is different, e.g. straight, angled or curved with large, oval-fenced, cross-scored mouth surfaces for a better grip. Some models may have replaceable, grooved rubber inserts in the jaw part. The instrument is made of stainless steel. It is a reusable instrument.

Placenta and Ovum Forceps; Uterine Elevating Forceps; Uterine Forceps: A surgical instrument used to grasp or manipulate the organ during an operation on the uterus. It is a reusable instrument.

Lung Grasping Forceps: A surgical instrument for atraumatically holding, manipulating or supporting the lung during a surgical procedure. It is typically scissor-like, self-holding with ring handles and is made of stainless steel. It's available in a variety of sizes and the gripping end is shaped like oval rings or triangular loops, which can be toothed to find better hold on the lung tissue. It is a reusable instrument.

Gall Stone Forceps; Kidney Stone Forceps; Lithotomy Forceps: A surgical instrument for gripping or manipulating urological stones to remove them. The instrument has two handles that need to be permanently compressed during use. For stone removal, other products may be used. It is a reusable instrument.

Tendon seizing instrument: A surgical instrument with a connected handle and two branches, usually serrated, for crossing, grasping, performing, holding or approaching sinews during surgery. It is a reusable instrument.

Bone Holding Clamps: A surgical instrument with stable branches and teeth for gripping or crushing bones. It is a reusable instrument.

Bone Holding Forceps; Cartilage Seizing Forceps: A surgical instrument for gripping and holding bones during a surgical procedure. It is usually sturdily constructed and has a scissor-like design that can be self-holding or non-latching, with ring handles or curved shaft handles. The instrument is made of stainless steel and is available in different sizes. The working end may consist of a wide variety of jaw variants, e.g., curved serrated jaw, straight, serrated jaw, or profiled key-shaped jaw with grooves for extra grip of the bones. This is a reusable instrument.

Hemostatic Forceps: A surgical device for applying temporary, atraumatic pressure to an artery except the artery arteriae pulmonales, ascending aorta, aortic arc, descending aorta, aortic bifurcation, coronary artery, common carotid artery, external carotid artery, internal carotid artery, arteriae cerebrales, truncus brachiocephalicus. For clamps used for a specific surgical purpose. It is a reusable instrument.

Vein Clamps: A surgical device for applying temporary, atraumatic pressure to a vein except the venae cordis, vena pulmonary, superior vena cava, and inferior vena cava. For clamps used for a specific surgical purpose. It is a reusable instrument.

Bone Cutting Forceps; Bone Rongeur Forceps: A surgical instrument with stable branches and teeth for gripping, cutting or crushing bones. It is a reusable instrument.

Cuticle- and Nail Instruments; Cuticle Nippers; Nail Nippers: An instrument for cutting fingernails and toenails. It can be bent to match the natural curvature of human nails. The toenail models are usually made more robust to cut thick toenails. The instrument consists of two movable branches, with holes for fingers and thumbs in the handle, and cuts by closing the leaves over the nail. (see also nail clippers).

Tissue- and Organ Grasping Forceps: A scissors-shaped, surgical instrument with ring handles whose working end is shaped as a ring, loop or ellipse to hold tissues or organs. It is a reusable instrument.

Catheter Introductory instrument: a metallic, flexible rod inserted into the lumen of a catheter or cannula to cause an increase in stiffness. This simplifies the manipulation and introduction into the body for various procedures. It is a reusable instrument.

Cartilage Crusher: A product by which autologous cartilage tissue is broken up. This material can be used to seal a perforation of the nasal septum. It is a reusable product.

Forceps, other: A surgical instrument for holding / gripping and / or squeezing tissue during a surgical procedure. It typically has two designs: 1) a self-retaining, scissors-like design with ring handles; or 2) a large tweezer-like design with straight branches interlocking at the proximal end. It is available in different sizes. The working end may be made in various designs, e.g. ring-shaped, triangular, with parallel grooving or toothing to give grip, or inwardly bent profiles for grasping. The industries are typically wide and thin. It is made of stainless steel. It is a reusable instrument.

Special pliers: An instrument with a hard metal inlay for removing FlexTIP tips. It has a scissor-like design with insulated ring handles. It is a re-usable instrument

3.2 Contraindications

There are no known contraindications.

4 Warning

!	Medical products are delivered in a non-sterile condition and must be cleaned, disinfected and sterilized prior to their initial use.
!	The use of faulty instruments is in principle forbidden and they have to go through the whole cleaning process before return.
!	Please take into consideration that through higher power a bigger damage of the tissue can result: f.e. on forceps: the power at the end of the jaw is higher than at the tip of the jaw
!	Please observe the additional information enclosed with the products.
!	Remove all protective sleeves and films prior to first using or preparation for use.
!	The safe combination of different products or of products with implants must be reviewed prior to clinical application by the user.
!	Avoid improper throwing or dropping of instruments
!	To prevent all contact corrosion, instruments with damaged surfaces must be separated immediately.
!	If used on patients with Creutzfeldt-Jakob disease or HIV infection, we are unable to accept any responsibility for reuse.

5 Handling

The type of treatment must be determined in each individual case by the surgeon in cooperation with the internist and the anaesthetist.

For operational use in various surgical disciplines must be done by appropriately trained and qualified personnel.

6 Preparation

The person in charge of preparatory treatment is responsible for ensuring that the treatment is duly carried out using the relevant equipment, materials and personnel in the treatment facility and so achieves the desired result. This necessitates validation and routine monitoring of the process used. We urge you to take note of the national regulations dealing with instrument preparation.

6.1 Reutilization restrictions

Frequent repeat preparatory treatment has minimal effects on the product. The end of the product life is normally determined by wear and damage due to use

6.2 Information on instrument preparation

- Use cleaning and/or disinfection agents with a pH-value within 9-10. Please observe manufacturer instructions regarding dosage, exposure time and renewal of solutions.
- Do not use hard brushes or coarse abrasive cleaners.
- Never leave instruments in cleaning or disinfection agents for longer than the specified time.

- Only used demineralized water for rinsing.
- Rinse and dry carefully through channels and pipes.
- Sensitive instruments must be cleaned in a storage or clamping fixture.
- Observe manufacturer instructions of cleaning – and sterilizing equipment.

6.3 Preparation at the place of use

Directly after using remove coarse dirt of the instruments and rinse out the working cannulas. Do not use fixing agents or hot water (> 40°C), as this results in residues becoming fixed and can affect the success of the subsequent cleaning operation

Dismantle and/or open instruments as far as possible. Within short time after use the instruments clean the instruments for reducing a drying of the residues. **Use cleaning and/or disinfection agents with a neutral pH-value <10.**

This enables an easier cleaning. If instruments come into contact with corroding medicines or cleaning agents, wash these up with water immediately after use.

6.4 Ultrasound bath (optional)

All instruments must be opened, dismantled and any cavities rinsed through.

Place instruments in the screen basket in such a way that overlaps and contact between instruments are avoided. Add cleaning

agent to the water and adjust the temperature of the solution in line with the cleaning agent manufacturer's instructions.

The cleaning in the ultrasound bath should be at 35-40 kHz, 5 minutes at least.

Subsequently rinse instruments including all cavities before cleaning and disinfection.

Medicine products which are possessing a bad ultrasound transmission, e.g. soft material are not usable for the ultrasound bath.

6.5 Manual cleaning

We recommend using lint-free soft cloths, plastic brushes and/or cleaning guns for manual cleaning. The cleaning agent should not have any fixing properties. The cleaning fluid must be changed daily.

Hints:

- If possible, apply the cleaning agent into cavities at pressure and carefully brush out.
- Be careful with the cleaning of barriers, slits, tight entries and joints etc.
- It has to be secured that instruments with joints will be cleaned in an opened and closed position.
- Use adequate brushes for the cleaning of cannulas, working channels and sack holes
- Flexible (complexe) instruments with unvisible surface and chuck have to be cleaned manually.
- We recommend for all push shafts, -and pipe shafts instruments and instruments whose surfaces are on top of each others during the cleaning (f.e.bone forceps and gouge forceps) a manually cleaning for an optimal cleaning result without residues

After that rinse the instruments with fully demineralized water extensively. Subsequently dip the instrument in disinfectant and use suitable cleaning instruments such as synthetic soft brushes.



Observe the cleaning agent manufacturer's instructions. At the end of the cleaning process, repeat the rinsing action with fully demineralized water and dry the instruments completely.

6.6 Mechanical cleaning



On the basis of international standards (EN ISO 15883) and national directives, only validated machine cleaning and disinfection methods may be used. For the mechanical cleaning we recommend a standard programme for surgical instruments, f.e. instruments from Miele.

Observe the following by loading:

- Place the dismantled/opened instruments securely in the tray.
- Instruments with openings and gaps have to be faced down with the opened side so that they can be cleaned and no water of the cleaning process is collecting inside them. If available use balanced devices for rinsing

- Place the instruments with joints in an opened position into the cleaning,-and disinfection machine.
- Do not overload trays, avoid creating any overlaps.

Preliminary rinsing (cold, if applicable fully demineralized water without additives) is followed by chemical.

The chemical cleaning should take place at 40°C -60°C for at least 5 minutes.

We recommend products with a **pH-value within 9-10**, e.g. MediClean forte from Dr.Weigert. The cleaning agents used should be selected depending on the material and properties of the instruments and in accordance to national regulations: If there is a high chloride concentration in the water, pitting and tension crack corrosion can occur on the instruments. The occurrence of this type of corrosion is minimized by using alkaline cleaning agents and demineralized water. By adding an acid-based neutralization medium, the rinsing off of alkaline cleaning agent residues is facilitated during the first intermediate rinsing process (warm or cold water).

In order to prevent the formation of deposits, it is advisable to use neutral cleaners where the water quality is unfavourable. After the second intermediate rinsing process, thermal disinfection takes place.

 The thermal disinfection should take place at temperatures of between 80 and 95°C, with an exposure time as outlined in EN ISO 15883.

After the finishing of the programme take the good out of the machine because corrosion can arise if the instrument remains in the machine.

6.7 Drying

Ensure adequate drying by the cleaning and disinfection device or using other suitable measures.

 The drying temperature should not exceed 95°C to avoid material-related ageing processes.

7 Sterilization

Prior to sterilization, products must undergo cleaning and disinfection, be rinsed off without residue using demineralized water and subsequently dried. HEBUmedical recommends using a validated steam sterilization process (e.g. sterilizer in compliance with EN 285 and validated in accordance with DIN EN ISO 17665-1).

 On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132° C)** with a **minimum dwell period of 4 minutes** and subsequent vacuum drying.

The vapour must be free of ingredients, recommended limiting values of feed water and vapour condensate are determined through EN 285.

Other sterilization processes are compatible but not validated from HEBUmedical.

When loading, observe the recommended total weight. After the sterilization, check the sterile product packaging for damage, and inspect the sterilization indicators

8 Maintenance, inspection

After cooling to room temperature, the instruments must be visually inspected for protein residues and other contamination, paying particular attention to cavities, blocks, inclusions, pipes, and other inaccessible areas. Instruments which are not free of residues must be returned for a complete retreatment process.

Carry out the functional check mentioned above. Instruments with stains, which are blunt, bent, no longer function or which are otherwise damaged must be segregated!

9 Packaging

Compliant packaging of products for sterilization in line with ISO 11607. Packaging used must be suitable for the instruments and protect them from microbiological contamination during storage. The seal must not be under tension. HEBUmedical recommends container for sterilized packaging.

10 Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from 5°C to 40°C. Protect from the effects of the sun's rays and artificial light. Warranty / Repair

11 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

12 Manufacturer and servicing address

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany
Tel. +49 7461 94 71 - 0
Fax +49 7461 94 71 - 22
eMail: service@HEBUmedical.de
Web: www.HEBUmedical.de

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**Our vast range of products covers over 10 000 various instruments.
They can be found in Internet or request our catalog.**

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